

COLLECTING THE EGGS EASILY AND SUCCESSFULLY WITH SENSE





If you're having trouble becoming pregnant, in vitro fertilisation (IVF) may be right for you. IVF is used to treat a range of fertility problems. An IVF cycle involves several different steps that will help you increase the chances of becoming pregnant. One of the first steps is collecting the eggs in order to begin the treatment.

The IVF journey and the egg collection step



Boosting the egg supply



Egg collection



Egg meets sperm, Fertilisation



Growing and selecting the best embryo




Transfer to the womb



The 'two week wait'

STARTING THE TREATMENT WITH YOU IN FOCUS

Egg collection may be an emotional and stressful experience for any woman going through IVF treatment. Your eggs will be collected through an ultrasound guided procedure where you may be sedated. A very thin needle is passed through the upper vaginal wall into the ovaries. Fluid, containing the eggs, is removed from each ovary under gentle suction.




Immediately after aspiration of the ovaries, the eggs are isolated from the fluid. The eggs are placed in a dish containing nutritious medium and then kept safely in an incubator. You may experience some cramps, feel a little sore or experience spotting (minor bleeding from the vagina).

DESIGNED FOR YOU

A successful egg collection needs to be fast, precise and retrieve the maximal amount of undamaged eggs without complications. A retrieval needle, called Sense, gives you just that.

Improved comfort

The unique design of Sense, with a thin tip and a larger needle body, provides ideal conditions for maximised control, precision, aspiration time and improved patient comfort during and after the procedure.¹ All at the same time with the same needle.



19%
Less pain

Less bleeding and pain

According to two studies, egg collection with the Sense needle gives significantly less bleeding compared to a conventional needle.^{1,2} Using Sense also reduces the overall pain experience. With Sense, patients experience 19 % lower pain compared to egg collections using a conventional needle.¹

IMPROVE YOUR EGG COLLECTION WITH SENSE



Egg collection with Sense gives you less overall pain, less bleeding, without prolonging the egg collection procedure or damaging the eggs.

Ask your clinician to find out more about Sense and how it can help you improve your chances for a successful IVF treatment.

REFERENCES

1. Wikland M et al. Pain experience during oocyte aspiration: Comparison between a newly designed reduced needle and a standard needle. Hum Reprod. 2010;25 (Suppl. 1) O-226
2. Nakagawa, et al., The effect of a newly designed needle on the pain and bleeding of patients during oocyte retrieval of a single follicle. J Reprod Infertil. 2015;16(4):207-211



WE SHARE YOUR DREAM

Vitrolife's vision is to fulfil the dream of having a baby. Our commitment to increase pregnancy rates has never been more dedicated. Together with equally devoted clinics, we are improving IVF success and fulfilling more couples' greatest dream. We are very proud to be a part of making this happen.

Find much more information about Sense on our website www.vitrolife.com

For professional use only.
Not intended for audiences in the USA.

Vitrolife 

Vitrolife · Box 9080 · SE-400 92 Göteborg · Sweden · Tel +46 31 721 80 00 · www.vitrolife.com